

Fact Sheet



Walking Aids



Walking Frame

A Lightweight Aluminium Folding Walking Frame with a wide base for maximum stability, ideal for home use. The lightweight aluminium walking frame features soft grip handles, one button folding function and height adjustable legs. This frame can be folded when not in use and is a real space saver. Ideal for travel use. It is available with or without wheels.

- * Lightweight Aluminium frame
- * 2 models (with or without wheels)
- * Easy push button fold
- * Non-slip rubber stoppers
- * Can be fitted with skis
- * Max weight user: 133kg



Walking Frame Accessories



Walker Ski

Product Code: SPW2T2

Net Bag

Product Code: VP179N

Walker Apron

Product Code: VP179T

Specifications—RM129F/RM129FW

Overall height.....	82-99	Unit weight.....	2.4/2.6
Overall depth	51	Weight capacity.....	133
Overall width.....	56/64		
Folded width	10		
Wheel size	4"		

Measurements in cm, weight in kg

Product Code: RM129F Unwheeled

Product Code: RM129FW Wheeled

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au