

Walking Aids



Walking Sticks - Folding

Rubber Handle Walking Stick



Product Code: Blue Lilac Green
Grey VP155SC VP155SA VP155SD
VP155SB

- * Height to handgrip 74-84cm
- * Hanging strap for easy storage
- * Easy fold
- * 5 height adjustments
- * Rubber foot to prevent slips
- * Max weight user: 100kg

Cane/Crutch Holder

Product Code: VP159C

- * Expandable
- * Suit most canes
- * Clips on table surface
- * Rubber base prevents slipping



Ergonomic Handle Walking Stick

- * Height to handgrip 71-94cm
- * Left or right hand versions

Product Code: Left Hand VP155
Right Hand VP156



Tri Ferrule

- * Size 19mm



Product Code: VS101

Quad Ferrule

- * Size 19mm



Product Code: VSPC2QT

Patterned Walking Stick



- * Height to handgrip

Black / Black Handle
Product Code: VP155F

Animal Print
Product Code: VP155JK

Blue Check
Product Code: VP155JP

Red Engraved
Product Code: VP155KB

Brown Engraved
Product Code: VP155KD

Blue Engraved
Product Code: VP155KE

Black Engraved
Product Code: VP155KG

Floral folding with
pattern handle
Product Code: VP155SG

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au