

Fact Sheet



Seating



Windsor Rise & Recline

The Windsor is a superior, luxurious and incredibly comfortable Theorem Rise and Recline Chair exclusive to Rehab and Mobility. With 4 motors providing an independent headrest, footrest extension and power lumbar support along with traditional rise and recline functions.

- * Modern, compact design
- * 4 individual functions
- * Independent head rest
- * Power lumbar support
- * Footrest extension
- * Blacklit easy to navigate hand control
- * USB connectivity on hand control
- * Side pocket
- * Lithium battery back up facility
- * 1 year warranty
- * Weight capacity P: 136kg S: 158kg



Specifications—Petite/Standard

Overall height..... 104/106
 Overall depth 98/100
 Overall width..... 84/89
 Seat depth 47/50
 Seat width 53/56

Seat height 47/49
 Armrest height 16/16
 Unit weight..... 68/71
 Weight capacity..... 136/158

Measurements in cm, weight in kg



Bark
Product Code:
 LRWP-589S-B Petite
 LRWS-589-B Standard

Fog
Product Code:
 LRWP-589S-F Petite
 LRWS-589-F Standard

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au