

Shower Chair detachable Back

The Affinity Shower Chair with detachable backrest offers the user the option of stool or chair. Durable anodized aluminium frame with plastic blown seat and backrest. Backrest has multiple positioning options, can be removed to use as a stool. Rubber armrests provide comfort and safety.

- * Lightweight aluminium frame
- * Removable backrest
- * Drain holes to reduce water build up
- * Height adjustable
- * Large non-slip stoppers
- * Weight capacity 113kg



Product Code: SCBAR

Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Unit Weight	Weight Capacity
75-87	43	51	41-44	41	41-54	15	3.5	113

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au