

Bathroom



Shower Stool w Padded Seat

The Affinity Shower Stool with soft PU seat offers the ultimate in comfort. Strong sturdy construction gives this chair an exceptional weight rating of 150kg. The soft removable PU seat pad makes cleaning a breeze. U shape cutout seat front assists with personal hygiene. Available in 3 seat widths and optional removable backrest.

- * Lightweight powder coated aluminium frame
- * Removable soft PU seat - easy cleaning
- * Drain holes to reduce slip and water build up
- * Height adjustable—5 choices at 2.5cm increments
- * Large non-slip stoppers
- * Available in 3 widths
- * Optional extra padded backrest
- * Weight capacity 150kg



Product Code: SC42.....42cm
SC47.....47cm
SC51.....51cm

Product Code: SC51BR

Specifications

	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat height	Armrest Height	Unit Weight	Weight Capacity
42	60-68	54	51	33	42	40-48	20	3.5	150
47	60-68	54	56	33	47	40-48	20	4.1	150
51	60-68	54	60	33	51	40-48	20	5	150

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au