

## Walking Aids *Bariatric*

### Bariatric Walking Frame

The Bariatric Walking Frame is a stable and helpful option for those needing additional support when walking. This frame is heavy duty, sturdy and height adjustable with a wide base. Foldable, making it easy for transportation and storage, bringing the sides in with two easy push buttons helps those with limited dexterity.



Product Code: RM804

- \* Sturdy heavy duty frame
- \* Height adjustable
- \* Easy push buttons to fold, transport and store
- \* Dual front and side brace for extra support
- \* Contoured foam handgrips for greater comfort
- \* Can be customized with walker skis
- \* Weight capacity 270kg

#### Walking Frame Accessories

VP179T Apron Bag

VP179N Net Bag

(details available [www.rehabandmobility.com.au](http://www.rehabandmobility.com.au))

#### Specifications

Overall Height	Overall Length	Overall Width	Folded Width	Unit Weight	Weight Capacity
80-105	52	72	11	5.7	270

Measurements in cm, Weight in kg

#### Your Local Agent

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [wholesale@rehabandmobility.com.au](mailto:wholesale@rehabandmobility.com.au)  
**NSW** Unit 2/144 Lake Entrance Road, Oak Flats 2529  
**Ph 1300 002 805** [nsw@rehabandmobility.com.au](mailto:nsw@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)