



Fact Sheet

Bathroom *Bariatric*

Bath Transfer Bench HDuty

The heavy duty durable aluminium Bath Transfer Bench is placed over the bath for those who have a shower over the bath situation. The wider seat allows the client to slide across from the outside of the bath. Two large suction feet for secure placement. Fully height adjustable will fit over most baths.

Optional extension legs available for external leg for higher bath tubs.

- * Durable aluminium frame
- * Moulded seat and removeable backrest
- * Heavy duty suction tips
- * Detachable armrest
- * Weight capacity RM425 280kg



Product Code: RM425

- * Optional extension legs available to raise bench into higher position

- **Complies with Australian Standards
AS/NZS 3973:2009**

Specifications

	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Tub Clearance Height	Armrest Height	Unit Weight	Weight Capacity
RM425	83-96	58	79	41.5	68	46-55	26	6.6	280

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au