



Bedroom



Bed Cradles

Bed Cradle is designed to keep blankets and weight away from legs and feet.

- * Keeps bed covers away from the feet and legs
- * Height adjustable for a comfortable fit to suit each individual
- * Can be fitted to the side or end of the bed as required
- * Suits any bed



Product Code: RM503



Product Code: VY440

Specifications RM503

Overall height 45-68
 Overall depth 72
 Overall width 36
 Unit weight 2.2
 Weight capacity 5

Specifications VY440

Overall height 45-68
 Overall depth 56 & 71
 Overall width 39
 Unit weight 2.5
 Weight capacity 5

Measurements in cm, weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
 www.rehabandmobility.com.au