

Coniston Swivel Bath Seat

The Coniston Swivel Bath Seat is a high quality rotating bath seat. Width adjustable frame and armrests, allows fitting to not only standard baths, but also extra wide baths (up to 79cm). The Coniston Bath Seat has a pivoting base that enables the seat to turn 90 degrees in either direction, it can be locked in position for maximum stability. The Coniston offers maximum stability and comfort making it an ideal alternative to many bath boards. No special fittings are needed for installation and can be removed in seconds.

- * Moulded plastic one-piece seat
- * All alloy frame
- * Width adjustable bars
- * Non slip stoppers
- * Width adjustable arms
- * Seat turn 360° in either direction
- * Four locking seat positions
- * Weight capacity 150kg



Product Code: VR130

Specifications

Overall Depth (mm) 610
 Seat Height (mm) 180
 Seat Width (mm) 440
 Overall Width (mm) 680-730
 Width between Arms (mm) 525-635
 Net weight 5.8kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
 Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au