



## Ellipse Heavy Duty Rollator

The Ellipse Heavy Duty Rollator is a high quality aluminium rollator perfect for use both indoors and outdoors. Featuring a wider seat width, height adjustable handles, a unique suspension system and a removable carry bag. This rollator offers smoother use over most terrains.

- \* Sturdy aluminium
- \* Height adjustable handles
- \* Unique suspension system
- \* 58cm seat width
- \* Perfect for indoor and outdoor use
- \* Removable storage bag
- \* Easy to fold, store and transport
- \* Locking clip to lock folded
- \* Weight capacity 205kg



### Specifications

Overall height.....	91-104	Folded width .....	23
Overall depth .....	69	Unit weight.....	10.4
Overall width.....	66	Weight capacity .....	205
Seat width .....	58		
Seat height .....	60		

*Measurements in cm, weight in kg*

**Product Code: 9346376000765**

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)