



## Wheelchair

### G4 Plus Wheelchair

The G4 Wheelchair is multi-adjustable allowing the chair to be tailored to the individual. Equipped with a number of important features along with a host of accessory choices, makes this a very desirable chair for everyday use. This model has an optional drum brake providing extra comfort for the passenger.

- \* Sturdy aluminium frame
- \* Height adjustable seat
- \* Height adjustable attendant assist brake
- \* Height and angle adjustable backrest
- \* Flip up, removable, height adjust, sliding, padded armrest
- \* Swing away removable leg rests
- \* Padded nylon upholstery
- \* Quick release 24" mag rear wheels with PU Tyres
- \* 7" front caster
- \* Optional drum brake
- \* Optional seat depth adjustment kit
- \* Width extendable foot plate
- \* Height adjustable Anti-tip bars
- \* Weight capacity 150kg



*Drum brake pictured*



Optional Elevating  
Leg Rest



Optional Stump  
Support



#### Specifications 46/51

Overall height .....	100-107	Seat height.....	51
Overall depth .....	106	Backrest height .....	48-53
Overall width .....	67/72	Armrest height.....	21-28
Seat depth.....	40-50	Unit weight .....	17
Seat width.....	46/51	Weight capacity .....	150

*Measurements in cm, weight*

**Product Code: RM144-46** 46cm Seat  
**Product Code: RM144-46D** Drum Brake  
**Product Code: RM144-51** 51cm Seat  
**Product Code: RM144-51D** Drum Brake

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)