

Scooter

Heartway Mirage Lux

Medium sized



Product Code: MS0020*

Enjoy the luxury of a supportive contoured seat, large chunky tyres and front and rear suspension to keep you comfortable as you visit the shops or explore your neighbourhood. A travel range of up to 30km gives you plenty of scope for outings.

The 360-degree swivel seat is height adjustable and the armrests, backrest and headrest can also be adjusted to suit your needs.

Large soft-touch buttons on the display panel assist with easy operation.



FEATURES

- Luxury comfort swivel seat
- Headlights, tail lights, deck lights, hazard lights, indicators
- Rear view mirrors
- Front carry basket
- Horn
- Speed control dial
- Armrests
- Built-in cupholders

SPECIFICATIONS

Motor Power	550w
Max Speed	10km/h
Operation Range	Up to 38km
Max Climbing Grade	6°
Turning Radius	1600mm
Ground Clearance	140mm
Max Loading Weight	136kg
Dimensions	1300x630x1220mm
Tyres	310 x 100mm pneumatic
Rear Caster (anti-tipper)	62x18mm
Weight with batteries	103.8kg
Weight without battery	72.2kg
Battery	12V/50ah x 2
Battery weight	31.6kg
Suspension	Front and rear