

Fact Sheet



Seating *Bariatric*

High Back Chair & Leg Rest

The High Back Orthopaedic chair offers maximum comfort and support with contoured backrest and padded armrests, perfect for post therapy recovery. Height adjustable along with seat depth adjustment has this chair suited to most users. Easy clean vinyl makes it ideal for both home or facility use. High back design provides increased postural support. Available in two seat widths & weight capacities



- * High back orthopaedic chair
- * Ideal for post therapy comfort
- * Strong sturdy construction
- * Powder coated steel frame,
- * Padded seat & armrests
- * Contoured padded backrest
- * Height adjustable
- * Dismantle for easy transportation
- * Easy clean Vinyl
- * Available in Grey
- * Weight Capacity 205kg
Bariatric 250kg



Specifications

Overall height..... 42-57
Overall depth..... 59
Overall width 37
Unit weight..... 4.3

Product Code: RM531LG

Specifications—Standard/Bariatric

Overall height..... 114-129	Seat height 47-63
Overall depth 65	Armrest height 15
Overall width..... 64/74	Unit weight..... 15/17
Seat depth 50-55	Weight capacity..... 205/250
Seat width 49/60	

Measurements in cm, weight in kg

Product Code: RM533 Standard
Product Code: RM533B Bariatric

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au