

Hythe Mobile Shower Chair

The Hythe Mobile Shower Chair incorporates the latest in seat design for improved safety and comfort. Height adjustable with 4 locking castors and features a pram handle as standard. The seat drainage outlets fully conform to current MHRA safety guidelines. Available with a variety of options for height adjustment and brakes.

- * Tubular steel frame with tough polymer coating
- * Pram handle for carer
- * Moulded polypropylene seat
- * 4 large lock in 5" castors
- * Easy to clean and maintain
- * Mobile shower chair ideal for assisting users in and out of the shower area
- * Can be used in wet environments and level access showers
- * Weight capacity 127kg



Product Code: VB560

Specifications

| Overall Height | Overall Length | Overall Width | Seat Depth | Seat Width | Seat Height | Armrest Height | Unit Weight | Weight Capacity |
|----------------|----------------|---------------|------------|------------|-------------|----------------|-------------|-----------------|
| 83-92 | 63 | 49 | 40 | 44 | 51-60 | 15 | 7 | 127 |

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au