

Fact Sheet



Wheelchair



Lightweight Aluminium Wheelchair

The Lightweight Aluminium Wheelchair is available in both self propel and transit options. With many of the features of other chairs the lightweight wheelchair offers ease of movement for the user and easy folding for trouble-free transportation.

- * Lightweight Self Propelled
- * Care assist hand brake
- * Folding backrest (for storage or transport)
- * Quick release wheels (for transport)
- * Flip up armrest for side transfer
- * 24" rear wheel, 6" front castor
- * Swing away and removable leg rests
- * Solid tyres (no maintenance required)
- * Weight capacity 110kg



Specifications AP/SP

Overall height	92/96	Seat height.....	50
Overall depth	103/107	Backrest height	43
Overall width	54/63	Armrest height.....	20
Seat depth.....	40	Unit weight	10/13
Seat width.....	46	Weight capacity	110

Measurements in cm, weight

Product Code: WH912-46 A/P

Product Code: WH910-46 S/P

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au