

Milton Aluminium Shower Chair

The Milton Stackable Shower Chair is a portable shower chair providing secure seat for those unable to stand in the shower. Incorporating the latest design in moulded shower seat for improved safety, drainage and comfort. Fixed armrests allow for safe transfer in and out of the shower.



- * Lightweight rust proof aluminium frame
- * Durable powder coat finish
- * Height adjustable
- * Moulded plastic seat with drain holes
- * Non-slip rubber tips
- * Weight capacity 150kg

• **Australian Standards**
AS/NZS 3973:2009

Product Code: VB499AL

Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Unit Weight	Weight Capacity
79-91	44	57	41	44	42-55	20	3.7	150

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au