

Aluminium Over Toilet Aid RM400

The RM lightweight aluminium over toilet aid with durable “S” bend side frame offers maximum comfort and space when seated. Supplied with privacy lid and large rubber suction ferrule to prevent slipping and provide maximum safety.



Product Code: RM400

- * Aluminium heavy duty frame
- * Included Splash guard
- * Privacy lid included
- * Easily removable seat
- * “S” bend side frame provides maximum space
- * Large rubber suction ferrule to prevent slipping
- * Weight capacity 160kg

Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Unit Weight	Weight Capacity
65-80	61	58	43	50	44-59	21	3.65	160

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au