

## Bathroom *Bariatric*

### Combo Shower Stool/Chair/OTA

The Combo Shower Stool/Chair/OTA offering open gap in front and rear seat allows access to body areas usually hard to reach. Can be used as a combination of an over toilet aid / commode or provide a secure safe option for those who require a seat in the shower. Suitable for home use, hospitals and nursing homes. Optional backrest available.



Product Code: RM371-49

- \* Lightweight rust proof aluminium frame
- \* Durable powder coat finish
- \* Height adjustable extension legs
- \* Suction rubber tips to reduce slip
- \* Soft polyurethane padded seat
- \* Includes bucket with lid and splash guard
- \* Very user friendly for both the carer and user
- \* Optional backrest
- \* Weight capacity 227kg



Overall Height	Backrest Width
86.5-98.5	49

Optional Backrest

Product Code: RM408-65BR

### Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Unit Weight	Weight Capacity
70-80.5	56	62.5	42	49	50-60	21	6.19	227

Measurements in cm, Weight in kg

### Your Local Agent

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [wholesale@rehabandmobility.com.au](mailto:wholesale@rehabandmobility.com.au)  
**NSW** Unit 2/144 Lake Entrance Road, Oak Flats 2529  
**Ph 1300 002 805** [nsw@rehabandmobility.com.au](mailto:nsw@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)