



# Fact Sheet

## Bathroom



### Pembury Corner Stool

The Pembury height adjustable corner shower stool, an excellent option when space is at a premium. Featuring fine leg adjustment maximises stability and safety when using the shower. A clip on/clip off seat with drainage holes allows for easy cleaning. Constructed from lightweight aluminium makes it easy to manage.



Product Code: VB636

- \* Ideal where space is at a premium
- \* Powder coated
- \* Height Adjustable
- \* Non-slip rubber tips
- \* Weight capacity 159kg

#### Specifications

Height (mm) .....	470-620
Seat Depth (mm) .....	265
Seat Width (mm) .....	530
Overall Depth.....	380
Maximum user weight (kg) .....	159

#### Your Local Agent

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085**    [wholesale@rehabandmobility.com.au](mailto:wholesale@rehabandmobility.com.au)  
**NSW** Unit 2/144 Lake Entrance Road, Oak Flats 2529  
**Ph 1300 002 805**    [nsw@rehabandmobility.com.au](mailto:nsw@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)