



Fact Sheet

Wheelchair Bariatric



Phonix Wheelchair

The European designed aluminium Phonix Wheelchair is an highly configurable self propel wheelchair. The flexibility and ease of adjustment in the seat width, seat height and seat depth offers a unique wheelchair for seat widths 50 - 65cm. Padded backrest and cushion seat offer the utmost in comfort.

- * High strength lightweight aluminium frame
- * Height adjustable attendant assist canes
- * Height adjustable, detachable armrest
- * Footrests angle adjustable
- * Full width foot plate option (as standard)
- * Swing away removable leg rests
- * 24" mag rear wheels with PU Tyres
- * Attendant assist hand brake
- * Anti-tip wheels (as standard)
- * Weight capacity 185kg



Product Code: PHONSK*

Optional extras:

- Elevating Leg rest
- Amputee support
- Tie down transport kit

Specifications

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest height	folded width	rear wheel	front castor	unit weight	weight capacity
92	108	70	42-54	50	42.5-52.5	14-24	36-46	17	32	24"	8"	25	185
92	108	73	42-54	53	42.5-52.5	14-24	36-46	17	32	24"	8"	25	185
92	108	76	42-54	56	42.5-52.5	14-24	36-46	17	32	24"	8"	25	185
92	108	79	42-54	59	42.5-52.5	14-24	36-46	17	32	24"	8"	25	185
92	108	82	42-54	62	42.5-52.5	14-24	36-46	17	32	24"	8"	25	185
92	108	85	42-54	65	42.5-52.5	14-24	36-46	17	32	24"	8"	25	185

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au