

Raised Toilet Seat with Armrests

The Raised Toilet Seat with Armrests is a convenient, ergonomically designed raised toilet seat. Designed to fit with majority of toilet bowl shapes. With tough PVC moulded seat made from light weight plastic for easy cleaning. Securely clamps into position on toilet pan and has 3 adjustable heights.



- * Ergonomically designed for comfort
- * Smooth aperture interior for easy clean
- * 3 height settings
- * Will fit most toilet suites
- * With lid
- * Soft padded flip up arm rest
- * Removable handle bars
- * Weight Capacity 136kg

Product Code: RM440

Specifications

Overall Depth (mm).....	470
Overall Width (mm)	575
Height Settings (mm).....	2"(50mm) 4"(100m) 6"(150mm)
Armrest Height (mm).....	200
Unit weight (kg)	4.85
Colour	White
Weight Capacity (kg)	136kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au