

Raised Toilet Seats

Raised Toilet Seats are a convenient, ergonomically designed raised toilet seat. Designed to fit the majority of toilet bowl shapes. Manufactured from high quality ABS recyclable plastic. The side mounted brackets of the raised toilet seat have rubber over-moulds and the front skirt has a rubberised pad for maximum grip and stability. Also available padded 4".



Product Code: VR224C

Product Code: VR224D



Product Code: VR224H



Product Code: VR224F

- * Ergonomically designed for comfort
- * Smooth aperture interior for easy clean
- * Large aperture and front cutaway
- * 3 heights 2" (50mm) 4" (100mm) 6" (150mm)
- * Will fit most toilet suite
- * With or without lid
- * Weight Capacity 190kg



Product Code: VR223P



- * Padded for extra comfort
- * Removable cover can be washed
- * Secure Velcro fastening



Product Code: VR224G

Specifications

Aperture Width (mm).....225
 Aperture Depth (mm).....270
 ColourWhite
 Seat Width (mm).....380
 Seat Depth (mm).....400

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
 Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au