



# Fact Sheet

## Wheelchair



## Ramps

### Lightweight Suitcase

- \* Suits majority of scooters and wheelchairs
- \* Provides access to vehicles and buildings
- \* Safety all-weather tread
- \* Easy to use and transport

#### Specifications

Folded width ..... 370mm  
 Channel width.....700mm  
 Weight capacity ..... 272kg  
 (WCR180)..... 300kg



Product Code: 910mm(3ft) **VA142C**  
 1520mm(5ft) **VA142K**  
 1830mm(6ft) **WCR180**

### Telescopic Channel

- \* Gritted surface for increased grip
- \* Lightweight

#### Specifications

Length collapsed (VA147F) .....744mm  
 (VA147S) ...1201mm  
 Channel width .....189mm  
 Weight capacity ..... 180kg



Product Code: 1220mm(4ft) **VA147F**  
 2130mm(7ft) **VA147S**

### Ramp Roll Up

Product Code: (3ft) **VA141**




- \* Easy access to houses and cars
- \* Portable & Compact
- \* Supplied with handy carry bag  
 (20x20x90cm)

#### Specifications

Length ..... 760mm  
 Width .....200mm  
 Height .....900mm  
 Weight capacity ...300kg

### Folding Suitcase

- \* Non slip surface
- \* Tri-fold for easy storage and transport



#### Specifications

Folded width .....360mm  
 Folded length.....630 & 915mm  
 Channel width .....740mm  
 Weight capacity ..... 272kg

Product Code: 1220mm(4ft) **VA143D**  
 1830mm(6ft) **VA143S**

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
 Ph 1300 368 085 [wholesale@rehabandmobility.com.au](mailto:wholesale@rehabandmobility.com.au)  
**NSW** Unit 2/144 Lake Entrance Road, Oak Flats 2529  
 Ph 1300 002 805 [nsw@rehabandmobility.com.au](mailto:nsw@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)