



Fact Sheet

Walking Aids

Rollator RM202 Modula

RM202 Rollator with height adjustable handles, available in 3 different seat heights provides a reliable rollator to suitable most users. Sturdy enough to stand up to the every day indoor and outdoor use, this would have to be one of the most popular rollators available.

- * Lightweight aluminium frame
- * Wider rear frame for maximum walking space
- * Available in 3 seat height options (X,Y,Z)
- * Ergonomic handles for comfort
- * Fabric padded seat, fabric bag
- * Weight user: 150kg



Specifications

Specifications X/Y/Z

Overall height....77-88/82-93/90-102	Seat height..... 50/55/59
Overall depth 69	Folded width..... 24
Overall width 64	Unit weight 7.6
Seat depth..... 33	Weight capacity..... 150
Seat width 44	

Measurements in cm, weight

Product Code:

RM202B-X Blue **RM202R-X** Red
RM202P-Y Purple **RM202R-Y** Red
RM202C-Z Champagne **RM202B-Z** Blue

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au