



## Walking Aids *Bariatric*

### Rollator RM204 & RM204Wide

RM204 Rollator constructed in tough heavy duty tubular aluminium frame, available in 2 widths. Featuring rear frame widened for maximum walking space, large capacity storage bag beneath a comfortable soft padded seat and a weight capacity of 200kg. Light enough to fold and put in the car.

- \* Lightweight aluminium
- \* Oval tubular frame construction
- \* Wider rear frame for maximum walking space
- \* Ergonomic handles for comfort
- \* Fabric padded seat, large fabric bag
- \* Weight user: 200kg



- **Australian Standards Tested**  
**ISO 11199-2:2005 (R2016)**

#### Specifications—Standard/Wide

Overall height .....	87-99	Seat height .....	54/55
Overall depth .....	73	Folded width .....	34
Overall width .....	62/71	Width between arms	45/52
Seat depth.....	36/37	Unit weight .....	9/9.5
Seat width .....	36/49	Weight capacity .....	200

*Measurements in cm, weight*

**Product Code: RM204R-Y Standard**  
**Product Code: RM204WB-Y Wide**

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)