

Scooter

Rothcare Boston

Pull-apart



Product Code: RCSKPL00113

This compact scooter is great for confined spaces, both indoor and out. It comfortably fits in the back of a small hatchback. The tiller folds down and the seat lifts off. Comfortable seat that is height adjustable to help you find the right driving position and swivels 360 degrees for easy mounting and dismounting.

FEATURES

- Pulls apart for transport or storage
- Height adjustable seat swivels 360° for easy mounting and dismounting
- Anti-tippers help minimise risk of tipping on corners or bumpy ground
- Removable basket
- Solid puncture proof tyres



SPECIFICATIONS

Motor Power	220w
Max Speed	8km/h
Operation Range	Up to 10-15km
Turning Radius	1200mm
Ground Clearance	80mm
Max Loading Weight	136kg
Dimensions	1020x610x900mm
Tyres	Solid puncture proof
Weight with batteries	48kg
Weight without battery	38kg
Battery	24V/12ah x 2
Battery weight	10kg
Charger	24V / 2A
Suspension	Front and rear suspension

