

### Skandia Over Toilet Aid

The Skandia combined toilet seat and toilet frame offers additional height for greater comfort and integral arms for easier transfer. Finished in durable hammertone powder coat helps prevent against scratches and scuffs, ideal for reissue (hire).



Product Code: VR157AL

- \* Lightweight rust proof aluminium frame
- \* Durable hammertone powder coat finish
- \* Height adjustable
- \* Supplied with splash
- \* White or grey seat and lid option
- \* Clip off seat for easy clean
- \* Non-slip rubber tips
- \* Pail & Lid available (sold separately)
- \* Weight capacity 150kg

• **Australian Standards**  
**AS/NZS 3973:2009**



\* Grey Pail w Lid for VR157AL  
Product Code: VS241

### Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Unit Weight	Weight Capacity
70-81	59	54	42	46	48-62	20	3.7	150

Measurements in cm, Weight in kg

### Your Local Agent

**QLD** 6-8 Imboon St, Deception Bay 4508  
Ph 1300 368 085 [wholesale@rehabandmobility.com.au](mailto:wholesale@rehabandmobility.com.au)  
**NSW** Unit 2/144 Lake Entrance Road, Oak Flats 2529  
Ph 1300 002 805 [nsw@rehabandmobility.com.au](mailto:nsw@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)