

SOLAX MALETA CARBON FIBRE

AUTOMATIC FOLDING MOBILITY SCOOTER



**LIGHTWEIGHT
CARBON FIBRE FRAME**



**REMOTE CONTROL
AUTOMATIC FOLDING**



**LITHIUM
BATTERIES**



SOLAX MALETA CARBON FIBRE

AUTOMATIC FOLDING MOBILITY SCOOTER

The Maleta Carbon Fibre is packed with new innovative modern design features. This model sees the addition of a lightweight carbon fibre frame, NFC keyless entry technology, a charging port on the dashboard, remote control for folding, headlights, ambient deck lights and ergonomic horn-shaped handlebars. It is safe, comfortable, easy to operate, convenient to be folded and carried, lightweight and steady. It can be folded and unfolded without any tools or effort.

SPECIFICATIONS

Weight With Batteries	20.8kg	Tyres	Solid - Puncture Proof
Battery Weight	2 x 0.85kg	Dimensions LxWxH (mm)	970x460x950
Weight Without Batteries	19.1kg	Folded Dimensions (mm)	280x460x650
Weight Capacity	125kg	Lights	Front, Rear and Ambient Deck
Adjustable Tiller Height	Yes	Drive Wheel	178mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance	25mm
Battery	24V/5.6ah Lithium x 2	Seat Dimensions	400mm x 330mm
Battery Range	Up To 15km	Seat Height From Ground	560mm
Maximum Speed	6km/h	Slope Ability	10°
Motor	120W	Turning Circle	1200mm
Number of Components	1 + Batteries	Suspension	No



6-8 Imboon Street
Deception Bay, QLD, 4508

1300 368 085
qld@rehabandmobility.com.au

Dealer Details