

# Fact Sheet



## Walking Aids



### Star 6" Rollator

The Star 6" rollator is a quick fit adjustable design meaning both the frame/seat and handles are height adjustable. This is great to get the perfect fit. Featuring a padded seat with flip up access to a secure privacy bag. This rollator is foldable and has a locking pin to keep it in its folded position for easier transport/storage.

- \* Height adjustable frame and handles
- \* Removable bag included
- \* Padded seat and back rest
- \* 6" castors
- \* Width between arms 48cm
- \* Max weight user: 150kg

Aust. STD Tested 493360



#### Specifications

Overall height .....	81-93	Seat height .....	49-57
Overall depth .....	65-73	Folded width .....	22
Overall width .....	62	Unit weight .....	7.3
Seat depth.....	35	Weight capacity .....	150
Seat width.....	37		

*Measurements in cm, weight*

**Product Code: V4261-FX**

**Your Local Agent**

QLD 6-8 Imboon St, Deception Bay 4508  
 Ph 1300 368 085 qld@rehabandmobility.com.au  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)