

Strood Aluminium Shower Stool

The Strood is a high quality shower stool designed for heavy duty use, with a maximum user weight of 159kg. Finished in durable hammertone powder coat helps prevent against scratches and scuffs.



Product Code: VB635AL

- * Lightweight rust proof aluminium frame
- * Durable hammertone powder coat finish
- * Height adjustable
- * Moulded plastic seat with drain holes
- * Clip off seat for easy clean
- * Non-slip rubber tips
- * Weight capacity 159kg

• **Australian Standards**
AS/NZS 3973:2009

Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Unit Weight	Weight Capacity
56-72	42	52	28	41	46-62	15	2.8	159

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au