

SWEETRICH TRACKER

LARGE SIZED MOBILITY SCOOTER



**CAPTAIN STYLE
SEAT & HEADREST**



**CANOPY AND
WINDSHIELD AVAILABLE**



**LARGE PNEUMATIC
TYRES & SUSPENSION**



SWEETRICH TRACKER

LARGE SIZED MOBILITY SCOOTER

The Sweetrich Tracker is a heavy-duty mobility scooter that features large 13-inch pneumatic tyres and offers a travel range of up to 60km in distance. The large captain style seat and rear suspension ensure a comfortable drive while catering for a 150kg weight capacity. Other handy features include a front rubber bumper, storage basket, cup holders, headrest and a digital dashboard.

SPECIFICATIONS

Weight With Batteries	128kg	Tyres	Pneumatic
Headlights	Yes	Dimensions LxWxH (mm)	1500x720x1200
Heaviest Part	107kg	Swivel Seat	Yes
Weight Capacity	150kg	Distance Between Armrests	520 - 620mm
Adjustable Tiller	Yes	Drive Wheel	330mm
Motor	800W	Ground Clearance	90mm
Battery	50ah Lead Acid	Seat Dimensions WxD	500mm x 480mm
Battery Range	Up To 60km	Seat Height From Ground	800mm
Maximum Speed	10km	Slope Ability	15°
Flip Up Armrests	Yes	Turning Circle	2200mm
Storage	Front Basket	Suspension	Front & Rear Suspension

Colours Subject to Availability



6-8 Imboon Street
Deception Bay, QLD, 4508
1300 368 085
qld@rehabandmobility.com.au