

# UNIVERSAL SWIVEL TRAY TABLE

Optional support handle making sitting and standing easy

Large bamboo 360degree swivel tray which also pivots out of the way when not in use

Height Adjustable

Adjustable base length

Large super grip rubber pads ensuring extra stability and protects flooring



## Available accessories



Handle



Organiser Compartment



Cup Holder



Side Utensil Tray

UNIVERSAL SWIVEL TRAY TABLE		
Height of Tray	61-81	cm
Height of optional handle	72-92	cm
Dimensions of optional handle	11.5x13	cm
Base Length Adjustment	65-90	cm
Weight	9	kg
Weight Limit of Tray	13.6	kg
Weight Limit of Handle	113.5	kg