

Walking Aids



Walking Frame

A Lightweight Aluminium Folding Walking Frame with a wide base for maximum stability, ideal for home use. The lightweight aluminium walking frame features soft grip handles, one button folding function and height adjustable legs. This frame can be folded when not in use and is a real space saver. Ideal for travel use. It is available with or without wheels.

- * Lightweight Aluminium frame
- * 2 models (with or without wheels)
- * Easy push button fold
- * Non-slip rubber stoppers
- * Can be fitted with skis
- * Max weight user: 133kg

Walking Frame Accessories



Specifications

Overall Height	Overall Length	Overall Width	Folded Width	Wheel Size	Unit Weight	Weight Capacity
82-99	51	56 64 with wheels	10	4"	2.4	133

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au