



Bed Rail & Cradle



Assembly

RM Bedroom products come ready assembled or may require minimal assembly.

Assembly and Adjustment of height (if applicable)

Step 1: Using your thumb, push down the Push-Pin and slide the frame up or down until the Push-Pin clicks back into the location hole.

(In some cases there maybe a push-pin on either side of the leg.)

Step 2: Ensure the Pin/s are completely through the hole, and test to ensure the two pieces do not come apart.

Step 3: Ensure both sides are adjusted to the same height.



Care & Maintenance

All of the product components are manufactured from high quality materials. They should be inspected regularly to ensure moving components (push-pins, locking screws) remain secure. If in doubt, do not use and immediately contact your supplier for support.

Your safety is of the utmost importance.

Warranty

All products carry a 12 month warranty from the date of purchase against manufacturing defects. Misuse or Damage is not covered. This warranty does not affect your statutory rights.

Your Local Dealer: